

## **X-GAL STAINING OF MOUSE EMBRYOS**

### **REAGENTS REQUIRED**

25% glutaraldehyde  
400mM potassium ferricyanide in ddH<sub>2</sub>O [K<sub>3</sub>(Fe(CN)<sub>6</sub>)], stored at -20 in the dark  
400mM potassium ferrocyanide in ddH<sub>2</sub>O [K<sub>4</sub>(Fe(CN)<sub>6</sub>)], stored at -20 in the dark  
20mg/ml X-Gal in DMSO or DMF, stored at -20 in the dark  
100mM EGTA, pH 8.0  
500mM MgCl<sub>2</sub>  
500mM Phosphate Buffer, pH 7.4  
10% Nonidet P-40 (NP-40)  
1% Sodium deoxycholate (Na-deoxy.)  
1M tris-HCl, pH 7.3  
Phosphate Buffered Saline

### **WORKING SOLUTIONS TO MAKE AND STORE AT RT**

β-Gal FIX: 0.1M PO<sub>4</sub>, 5mM EGTA, 2mM MgCl<sub>2</sub>

100mL 0.5M PO<sub>4</sub> buffer  
25mL 100mM EGTA  
2mL 0.5M MgCl<sub>2</sub>  
373mL ddH<sub>2</sub>O

β-Gal WASH: 0.1M PO<sub>4</sub>, 2mM MgCl<sub>2</sub>, 0.01% Na-deoxy., 0.02% NP-40

100mL 0.5M PO<sub>4</sub> buffer  
2mL 0.5M MgCl<sub>2</sub>  
1mL 10% NP-40  
5mL Na-deoxy.  
392mL ddH<sub>2</sub>O

β-Gal STAIN: 0.1M PO<sub>4</sub>, 2mM MgCl<sub>2</sub>, 20mM tris-HCl pH 7.3, 0.01% Na-deoxy, 0.02% NP-40 (+5mM K<sub>3</sub>(Fe(CN)<sub>6</sub>), 5mM K<sub>4</sub>(Fe(CN)<sub>6</sub>), on 1mg/ml X-gal)

100mL 0.5M PO<sub>4</sub>  
2mL 0.5M MgCl<sub>2</sub>  
10mL 1M tris-HCl, pH 7.3  
5mL 1% Na-deoxy.  
1mL 10% NP-40  
344.5mL ddH<sub>2</sub>O

Make up solutions minus last three components which will be added immediately before use

6.25mL 400mM K<sub>3</sub>(Fe(CN)<sub>6</sub>)  
6.25mL 400mM K<sub>4</sub>(Fe(CN)<sub>6</sub>)  
25mL 20mg/ml X-gal

Protocol:

- 1 Dissect out embryos into a dish of cold PBS
- 2 Wash embryos once in cold PBS on ice (~5 mins without rocking)
- 3 **Fix the embryos**  
Add 400 $\mu$ L glutaraldehyde to 50mL cold FIX solution. Aliquot 1.2mL cold fixative into 2mL round bottomed microtubes containing one embryo each.  
Incubate on ice for:
  - E6 to E8 – 10 to 15 min
  - E9 to E11 – 15 to 20 min
  - E12 to E13.5 – 20 to 30 min
  - \*\* for E13.5, use 25 min
- 4 **Wash fixed embryos.**  
Transfer embryos to pre-labelled 12 well dish. Incubate 3x 5min in 2mL of  $\beta$ -gal WASH solution at room temperature with rocking.
- 5 **Stain embryos**  
Pre-warm stain working solution to 37°C. Just before adding stain to wells, add the appropriate volume of K3(Fe(CN)6), K3(Fe(CN)6), and X-gal.  
Add 2mL stain per well.
  - Stain for 1 hour to overnight (or longer) at 30-37°C. Avoid incubators with added CO<sub>2</sub>.
  - \*\* Tris-HCl is only required in stain solution if staining overnight
  - \*\* If incubating for longer than a few hours, be sure to wrap plate in parafilm
  - \*\* Both specific staining and background staining increase with increased incubation time.
  - \*\* When looking for activity in melanocytes, embryos are stained at 37°C for 48 hours
- 6 **Transfer embryos to a fresh 12 well dish and wash stained embryos 3x in PBS at RT with rocking for (2 minutes for first wash, 5 minutes for 2<sup>nd</sup> and 3<sup>rd</sup> wash)**  
Pictures can be taken at this point or fixed and photographed later
- 7 **Fix stain with 4% formaldehyde in PBS at RT for at least 1 hour (with rocking)**  
Embryos can be stored in this solution
- 8 **Wash embryos in PBS (2x 5 minutes at room temperature with rocking) and photograph**
- 9 **Clearing**  
Incubate embryos in 50% glycerol, 50% PBS, then  
80% glycerol, 20% PBS
  - \*\* incubate in glycerol:PBS mixture until embryo has gone from floating on top of solution to being submerged and laying on bottom of tube or well
  - \*\* can store in the last glycerol:PBS solution
- 10 **Genotyping embryos**  
Cut off a small portion of the bottom of the embryo (ex. Tail + 1 leg) and digest using proteinase K.
  - \*\* Digestion is more successful after embryos have been cleared vs. when sitting in formaldehyde fix solution