

REVERSE TRANSCRIPTION – PCR OF MOUSE TISSUE RNA

REAGENTS REQUIRED

TRIzol (Invitrogen, Cat.15596-026)
Chloroform
Isopropanol
75% Ethanol (in DEPC-treated water)
DEPC-treated water or RNase free water
DNase I Amplification grade (Invitrogen, Cat. 18068-015)
Promega GoScript Reverse Transcription System (Fischer Cat. PRA5000)
Regular PCR Reagents (5x GoTaq buffer, dNTPs, spermadine, primers, Taq, ddH₂O)

- 1 Homogenize mouse tissue
Grind up one whole ear or 1cm² piece of shaved flank under liquid nitrogen using a mortar and pestle previously chilled with liquid nitrogen. Transfer these skin grounds or 0.5cm³ sized other tissue into a 14mL disposable falcon tube and add 1mL TRIzol. Homogenize tissue using the handheld homogenizer at approximately half speed until tissue has been fully homogenized. Clean/rinse homogenizer between samples by placing end of homogenizer at full speed into three different solutions sequentially: 2 flasks with ddH₂O and one flask with 70% EtOH. Expel EtOH on outside of homogenizer by running homogenizer at full speed in the air/within the hood.
- 2 Extract RNA using TRIzol reagent following manufacturers instructions.
Resuspend RNA pellet in 30-50µL depending on pellet size. Final RNA concentration should ideally be between 150-1000ng/µL. If the concentration exceeds 1000ng/µL, dilute the sample to below this in order to effectively treat the sample with DNase.
- 3 Treat RNA samples with DNase
Treat 1µg RNA in a 10µL volume. Treat RNA following manufacturers instructions paying close attention to the length of each incubation step.
- 4 Reverse transcribe the RNA (make cDNA) using all reagents from Promega GoScript Reverse Transcription System
*NOTE – The following protocol is the one I used. Refer to the manual for troubleshooting options.
Only complete +/- RT for important samples (there is not enough buffer to complete this for all samples)

Part 1: RNA/primer pair

Add 4µL of DNase treated RNA to 1µL (0.5µg) random primers in 8 well PCR strip tubes
-random primers work better than the polydT primers
-use maximal amount of RNA (ie. 4µL) because it is at a lower concentration than can be used
-same volume of all samples because concentration already normalized by the amount of RNA treated with the DNase

→ in a PCR machine, incubate RNA/primer pair at 70°C for 5min, drop temp down to 4°C, then chill on ice

Part 2: RT mix

- use all solutions that come in the kit
- make a master mix using the following volumes of each reagent per RNA sample:

7.3 μ L nuclease-free H₂O
4 μ L GoScript 5X buffer
1.2 μ L MgCl₂ (25mM)
1 μ L dNTP
0.5 μ L rRNasin RNase inhibitor
1 μ L RT
15 μ L

→add 15 μ L RT mix to each well (containing primed RNA), spin down, and start the RT cycle
25°C – 5min
42°C – 60min
70°C – 15min
4°C – 10sec

5 PCR of target mRNA

Use 1 μ L of cDNA per PCR reaction in a 25 μ L volume.

Use 0.75 μ L of Taq made in house or 0.2 μ L (1 unit) of Promega GoTaq

17.4 μ L ddH₂O
5 μ L 5x GoTaq buffer
0.3 μ L spermadine
0.25 μ L dNTP (each dNTP at a concentration of 10 μ M)
0.3 μ L primers (each primer at a concentration of 10 μ M)
0.75 μ L Taq (lab-prepared)
24 μ L

To increase cDNA detection level, a nested-PCR can be employed. In this case, take into account the amount of buffer already contained within the PCR product being added as a template for the second PCR. Normally, I add 2 μ L of product from the first PCR as a template for the second PCR. In this case, the PCR mix per reaction would resemble the following:

16.8 μ L ddH₂O
4.6 μ L 5x GoTaq buffer
0.3 μ L spermadine
0.25 μ L dNTP (each dNTP at a concentration of 10 μ M)
0.3 μ L primers (each primer at a concentration of 10 μ M)
0.75 μ L Taq (lab-prepared)
23 μ L

6 Visualize PCR products by gel electrophoresis

Load PCR products onto a 1.8-2% agarose gel. The 20 well comb will provide sufficient sample separation.